



# KINSA FOR SCHOOLS

## TRAINING MATERIAL FOR PARENTS





# Initial Set-Up



# Download the Kinsa App

You can download the Kinsa app on:

- the [App Store for iOS](#)
- or [Google Play for Android](#)



**Kinsa** 12+

Health & Symptoms Guide  
Transform Health Inc.

#45 in Medical

★★★★★ 4.8, 10K Ratings

Free



## Kinsa for Wireless Smart Thermometers

**Kinsa Inc.** Health & Fitness

**E** Everyone


 This app is compatible with your device.



# Create your Account on the Kinsa App

WELCOME TO KINSA

## Helping families feel better faster



**Know what to do next**  
Care recommendations specific to your fever, age and symptoms.

• • •

**SIGN UP**

**LOG IN**

←

YOUR PRIVACY

## Privacy Policy

**Privacy & Disclaimer**

- We never use or share your personal information without your permission.
- We are not a substitute for your doctor, and our guidance is not appropriate for everyone.

I agree to the following:

- [Terms and Conditions](#)
- [Privacy Policy](#)
- [Medical Disclaimer](#)

**AGREE**

←

LET'S GET STARTED

## Create your account

email

password

Password must be at least 8 characters.

Enable fingerprint sign in


**SIGN UP**

Have an account? [Log in here.](#)

GREAT!

## Create your profile

Track your temperatures and receive age-specific fever guidance.



Name  
Mia

Birthday  
Jun 17, 1995



Gender  
Female

**SAVE**



# Add your Family Members (Optional)




- You have the option to add your family members when creating your account
- Your family members won't have to enroll to the Kinsa For Schools program (this is only for students and staff)
- You can always add, edit, or remove family members profiles in Settings (the gear icon on the top right of the home screen)





# Pair your Thermometer to the Kinsa App

- After creating your account, you will see this “Let’s get started!” card on the top of your feed
- Tap on “SET UP KINSA” and turn your thermometer on
- Your thermometer will be connected to the Kinsa app via Bluetooth: make sure that your Bluetooth is turned on as well
- Note that you can also pair your thermometer in Settings (gear icon on the top right of the main screen)






# SCHOOL SETUP



# Step 1: Join Your School Group

1. Tap Settings gear in the top right corner, then click 'Join my school group'
2. Search your school's name
3. After selecting your school, select the family members that you would like to add to the school. Specify their School Role and Grade and tap SAVE







# Daily Self-Check



# Daily Self-Check

- You will be prompted to complete a self-check each day
- As part of the daily self-check, you will be asked to:
  - Take your temperature, using your Kinsa thermometer
  - Answer a question around COVID-19 exposure
  - Select any symptoms you may have (if any)

**Would you like to take {Name}'s temperature?**

The CDC recommends checking your child's temperature each day before heading to school.


**TAKE TEMPERATURE**

**NOT NOW**



# If you do not pass the screen


- If you do not pass the daily screening you will receive a stay home from school card in the app.





# If you pass the screen and mark as feeling 'well'



- If you pass the daily screening you will receive a fun 'high-five' card in the app, signaling you're feeling well!





# If you'd like to share your information with others (for example your school)

After taking your temperature and entering your symptoms click on 'See Full Timeline'.



In the full timeline view a click the share button in the top right corner and select the share history.



Thank you for participating in  
Kinsa For Schools !